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Email : [CustomerService@BattleCreekEquipment.com](mailto:CustomerService@BattleCreekEquipment.com)

Adresse Postale : Battle Creek Equipment Company  
702 S. Reed Road  
Fremont, IN 46737

 **BATTLE CREEK**  
Equipment  
Enhancing Comfort • Enhancing Treatment • Enhancing Lives™

702 S. Reed Road / Fremont, IN 46737  
800-253-0854 [www.BattleCreekEquipment.com](http://www.BattleCreekEquipment.com)

© Battle Creek Equipment Company 2012 Rev C  
PN 3029

# Thermophore® Automatic Moist Heat Pack™

## Instructions

Please read and follow all instructions in this booklet,  
as well as those on the box and pad

 **BATTLE CREEK**  
Equipment  
Enhancing Comfort • Enhancing Treatment • Enhancing Lives™

702 S. Reed Road / Fremont, IN 46737  
800-253-0854 [www.BattleCreekEquipment.com](http://www.BattleCreekEquipment.com)

## **IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY**

### **DANGER**

To reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions:

- 1. READ AND FOLLOW ALL INSTRUCTIONS IN THIS BOOKLET, ON THE BOX AND ON THE PAD.**
- 2. DO NOT USE WHILE SLEEPING.**
- 3. DO NOT USE ON AN INFANT.**
- 4. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, PARALYZED PERSON, SLEEPING OR UNCONSCIOUS PERSON, OR A PERSON WITH POOR BLOOD CIRCULATION UNLESS CAREFULLY ATTENDED.**
- 5. DO NOT USE ON AREAS OF INSENSITIVE SKIN.**
- 6. THE MOMENTARY SWITCH WAS DESIGNED FOR YOUR SAFETY. USE ONLY AS DIRECTED ON PAGE 2.**
- 7. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING, CHECK SKIN UNDER PAD FREQUENTLY.**
- 8. NEVER USE PAD WITHOUT COVER IN PLACE.**
- 9. DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.**
- 10. DO NOT SIT ON OR CRUSH PAD-AVOID SHARP FOLDS.**
- 11. NEVER PULL THIS PAD BY THE SUPPLY CORD. DO NOT USE THE CORD AS A HANDLE. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.**
- 12. CAREFULLY EXAMINE INNER COVER BEFORE EACH USE. DISCARD THE PAD IF INNER COVERING SHOWS ANY SIGN OF DETERIORATION.**
- 13. THIS UNIT IS NOT DESIGNED FOR USE IN CONJUNCTION WITH ANY TYPE OF MASSAGE APPARATUS. SUCH USE MAY CAUSE HAZARDOUS DAMAGE TO INTERNAL COMPONENTS AND WILL VOID THE WARRANTY.**
- 14. SAVE THESE INSTRUCTIONS.**

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**Battle Creek Equipment Company**  
702 S. Reed Road  
Fremont, IN 46737

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#### **USING YOUR THERMOPHORE® MOIST HEAT PACK**

1. Slide the heating unit into the special cotton-blend fleece outer cover. **DO NOT** use the heating unit without the fleece cover, and do not add additional layers of insulating cover.

**DO NOT** add any water to any part of the Thermophore® Heat Pack. The outer cover is designed to absorb moisture from the atmosphere between use, and release the moisture when the heat is activated. Your skin should feel slightly damp from the water vapor condensed during the treatment, but will not be dripping wet. If you feel your cover is no longer producing moisture, you may simply need to wash it, or to let it rest in an area of high humidity for a few hours.

2. It is very important to keep the heating unit flat inside the cover. Folded or overlapping areas of the heating unit could intensify the heat. This could damage the unit, burn the user or cause a fire. Feel inside the cover with your hand to make sure the heating unit is flat and wrinkle-free. Before each use, check the heating unit to be sure it has remained flat inside the cover.
3. Insert the plug into any convenient 110-volt electric outlet. This heating pad has a polarized plug (one blade is wider than the other) as a safety feature. This plug will fit into a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the orientation of the plug in the outlet. If it still does not fit, contact a qualified electrician. Do not attempt to defeat this safety feature.
4. When you are ready to begin your therapy session, position the Thermophore® Heat Pack over the area you wish to treat. Direct contact with bare skin produces the most benefit. Make sure the unit is lying smoothly on your body so it can conform to your body's curves. **Do not lie or sit** on the Thermophore® Pack.  
The Velcro® strap on the Petite model will allow secure positioning.

#### **If you have the original, hand-held switch:**

5. Press the lever on the hand control and hold it down. **WARNING:** This switch was designed for your safety - **USE ONLY AS DIRECTED.** The switch is spring-loaded so the Thermophore® Pack is "off" if the switch isn't held "on". You control the heat by holding the switch in your hand. Please do not switch it "on" any other way.

When the Thermophore® Pack gets as hot as you desire, release the switch until you begin to notice that the unit feels cool. Repeat this heating and cooling cycle as many times as you like.

### **If you have the One-Touch hands-free timer switch:**

The One-Touch™ switch is a user-friendly electronic timer and temperature control. The small red light at the top of the switch will remain lit whenever the heating unit is plugged in properly. When you are ready for a surge of heat, just press the larger red “Power” button in the center of the switch display.

All four green bars should light up in a blinking cascade rotation to indicate a total of 20 minutes, divided into in four 5-minute time segments. One light-bar will become a solid green light each time the “Power” or “Resume” button has been activated for a total of five minutes, while the remaining bars continue to blink. This will help you monitor the length of your session, but you can end the session at any time, or extend it longer if you prefer, with the “Power” button.

As soon as the MaxHEAT™ therapy pad becomes as hot as is comfortable for you, just press the yellow “Pause” button, located directly below the “Power” button. This feature allows you to customize the heat level by allowing the heating pad to cool a bit while the timer function pauses, your body adjusts to the heat level and your tissues begin to relax. The “Pause” button will blink yellow, and the cascading green lights on the timer will pause. Whenever you feel that more heat is needed, push the yellow button a second time to “Resume” heating. The MaxHEAT™ pad is designed to get very hot, very quickly. Please stay alert, monitor the heat level, and do not let the heat level become uncomfortably hot.

If the “Pause” feature is left in place for 20 minutes, or the heating function has been activated for a total of 20 minutes by using either the “Power” or “Resume” function, the unit will automatically shut off as a safety feature. If you feel that you need extended heating at that point, just press the red “Power” button again to activate a new 20-minute timer cycle.

Because different people have different needs, the time for the entire session may vary from person to person. We recommend that the session not last over 30 minutes, with a minimum cooling period of 60 minutes between heat therapy sessions. Please consult your physician if you require a more intense heat session. Should unusual symptoms appear, discontinue use of the product and consult your physician.

When finished with your therapy session, remove the plug from the outlet. Store the Thermophore® Heat Pack in the original box or lay it out flat. Loop cord loosely when storing. Tight wrapping may damage cord and internal parts.

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If you live in an area with low atmospheric moisture, you may want to hang the fleece cover in an area of high humidity so that it can re-absorb moisture to be ready for the next session.

## HOW DOES THE THERMOPHORE AUTOMATIC MOIST HEAT PACK PRODUCE MOISTURE?

Most people are amazed to find that they can have an effective, penetrating moist heat treatment without adding water to their unit. The moisture is produced by drawing humidity from the air and retaining it in the outer fleece cover. When the Thermophore<sup>®</sup> Pack is turned on, the rapid heat rise forces moisture out of the fleece cover and onto the user's skin. A lot of moisture is not necessary—only a light amount is needed to transmit the intense heat more efficiently.

You can confirm that your unit is producing moisture by laying it on a hard, flat, heat-resistant surface (like a countertop or a piece of glass). Activate the switch, and allow the Thermophore<sup>®</sup> Pack to heat for about 5 minutes. Lift the unit and check underneath. You should be able to feel and see tiny beads of moisture on the fleece cover and probably on the surface, too.

If you feel that your Thermophore<sup>®</sup> Automatic Moist Heat Pack is not producing enough moisture, it may be because it has been used for an extended period and the fabric has dried out. You may also have extremely dry air because of your climate, air conditioning or heating system. **Do NOT add water.** Instead, try these suggestions:

1. Wash the outer fleece cover and either hang it to dry or use the medium setting on your dryer, so that it does not become excessively dry.
2. Avoid storing your Thermophore<sup>®</sup> Pack in a hot, dry location.
3. Allow time between heat therapy sessions for the cover to regenerate and absorb moisture.
4. Hang the cover near your shower or laundry, where the air contains more moisture. Caution: Do *not* put the inner heating unit of the Thermophore<sup>®</sup> Automatic Moist Heat Pack in or near water.

**WARNING** - USE YOUR THERMOPHORE<sup>®</sup> AUTOMATIC MOIST HEAT PACK<sup>TM</sup> CAREFULLY. FOR IT TO WORK PROPERLY, IT MUST BE CAPABLE OF BECOMING VERY HOT. PLEASE STAY ALERT. DO NOT USE OVER INSENSITIVE SKIN OR IN THE PRESENCE OF POOR CIRCULATION. PERSONS WITH DIABETES SHOULD BE ESPECIALLY CAREFUL WHEN APPLYING HEAT. THOSE WHO ARE INCAPACITATED SHOULD BE CLOSELY SUPERVISED IN THEIR USE OF THE THERMOPHORE<sup>®</sup> PACK AND SHOULD

NOT BE LEFT UNATTENDED. DO NOT USE A THERMOPHORE<sup>®</sup> AUTOMATIC MOIST HEAT PACK AS A BEDWARMER OR FOOTWARMER. IT GETS TOO HOT TO BE USED FOR THIS PURPOSE. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS. DO NOT USE YOUR THERMOPHORE<sup>®</sup> AUTOMATIC MOIST HEAT PACK ON ANIMALS OR SMALL CHILDREN.

### FULL TWO-YEAR WARRANTY

During the first two years after purchase, we will replace this product free of charge to the owner if it has a manufacturing defect. Please note that this warranty does not cover incidental or consequential damages or costs connected with any such defect. While this warranty gives you specific legal rights, you may also have other rights which vary from state to state.

Please be sure to save your sales receipt showing the date of this purchase. If no receipt is available, we will use the serial number on the unit to determine warranty eligibility.

### SERVICE CHECKS

If your Thermophore<sup>®</sup> Heat Pack fails to heat at any time, make these checks:

1. Make sure the plug fits properly and firmly into the electric outlet, and that the outlet is functioning properly. Try plugging the unit into another outlet which is on a different electrical circuit. Check your fuses or circuit breakers, if necessary.
2. Inspect the electrical cord. Make sure there is no break in the wires of the cord and no looseness of the cord next to the plug, switch, or entrance to the Thermophore<sup>®</sup> Pack. If the electrical cord has been damaged, it should only be serviced by the Battle Creek Service Department.
3. Check the heating unit inside the fleece cover. The duck cover of the heating unit may eventually become somewhat discolored, especially near the wires. This is to be expected. Excessive scorching or isolated darker areas caused by use of the unit in a folded or insulated condition, however, may indicate damage to internal components. This damage is not a manufacturing defect and will not be covered by our warranty. If this has occurred, the unit should be replaced.

If warranty service is required, or you have questions about the use of your Thermophore<sup>®</sup> unit, please contact our company for instructions.

**Battle Creek Equipment Company**  
702 S. Reed Road  
Fremont, IN 46737  
800-253-0854 or [www.BattleCreekEquipment.com](http://www.BattleCreekEquipment.com)

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