# Thermophore® MaxHEAT™

Heat Pack for Deep-Heat™ Therapy

## Warranty

During the first two years after purchase, we will replace this product free of charge to the original owner if it has a manufacturing defect. Please be sure to have your sales receipt showing the date of purchase. If no receipt is available, we will use the serial number stamped on the unit to determine warranty eligibility. Please contact us for assistance at 800-253-0854 or at CustomerService@BattleCreekEquipment.com. In order to obtain warranty service, you will be asked to send the unit to us for evaluation.

\* Please note this warranty does not cover incidental or consequential damages or costs connected with any such defect. While this warranty gives you specific legal rights, you may also have other rights which vary from state to state.

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## Instruction Manual Manuel d'instructions

### Constant-Comfort™ Switch

Three temperature levels let you select the temperature you prefer

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Meets or Exceeds  
Standards for UL130  
Certified to CAN/  
CSA C 22.2

Enhancing Comfort, Enhancing Treatment, Enhancing Lives!™

702 South Reed Road  
Fremont, IN 46737  
800-253-0854

www.BattleCreekEquipment.com

## IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

### DANGER

To reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions:

READ AND FOLLOW ALL INSTRUCTIONS IN THIS BOOKLET, ON THE PACKAGING AND ON THE PAD.

DO NOT USE WHILE SLEEPING OR IN BED.

DO NOT USE ON AN INFANT.

THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, PARALYZED PERSON, SLEEPING OR UNCONSCIOUS PERSON, OR A PERSON WITH POOR BLOOD CIRCULATION UNLESS CAREFULLY ATTENDED.

DO NOT USE ON AREAS OF INSENSITIVE SKIN.

BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING - CHECK SKIN UNDER PAD FREQUENTLY.

NEVER USE PAD WITHOUT COVER IN PLACE.

DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.

DO NOT SIT ON OR CRUSH PAD-AVOID SHARP FOLDS.

NEVER PULL THIS PAD BY THE CORD. DO NOT USE THE CORD AS A HANDLE. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.

CAREFULLY EXAMINE INNER COVER BEFORE EACH USE. DISCARD THE PAD IF INNER COVERING SHOWS ANY SIGN OF DETEIORATION.

THIS UNIT IS NOT DESIGNED FOR USE IN CONJUNCTION WITH ANY TYPE OF MASSAGE APPARATUS. SUCH USE MAY CAUSE DAMAGE TO INTERNAL COMPONENTS AND WILL VOID THE WARRANTY.

SAVE THESE INSTRUCTIONS.

## Using Your Heating Pad

- Insert the plug into any convenient 110-volt outlet. The pack has a polarized plug as a safety feature so the plug will only fit into a polarized outlet one way.
- Slide the control button to the heat level you prefer (High, Medium, or Low.) The small light at the top of the switch will glow as a reminder that the unit is On. Slide to position O to turn Off at any time. The switch will automatically turn off after approximately 25 minutes of continuous use. You can adjust the temperature at any time during your therapy session.
- We recommend treatments of 30 minutes or less with a minimum wait between sessions of 60 minutes. Please consult a physician if you require more prolonged heat therapy or if you notice unusual symptoms.

## Cleaning Instructions

- Do not soak the heating unit. Wipe clean only.
- Machine wash the Fleece Cover in cold water on the gentle cycle. Hang dry or dry flat.

## Storage Instructions

- Store your heating pad flat or in its original packaging.
- Do not wrap cord around pad for storage.

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